


Grand Rapids, MI 49512 • 616.656.2900

SS2712 - DUAL SINK SS2713 - TRIPLE SINK





#### INTRODUCTION

**SKYTRON Clean Up Stations** are used for the pre-cleaning, inspection and sorting during the decontamination process of surgical instruments.

### **TECHNICAL DATA**

**SKYTRON Clean Up Stations** are available in two models, a two bay sink and a triple bay sink. SKYTRON Clean Up Stations have sound deadening material added to the underside of the countertop, sink bottom and backs. Sink basins are 18" x 24" x 12"D (457 x 610 x 304mm), backsplash is 8" high with 45° edge. To each side of the basins there are pitched drain board to alleviate the pooling of water as well as a rolled edge on the sides and front to retain any water within the clean up station perimeter.

#### STANDARD ACCESSORIES

Each SKYTRON Clean Up Station shall include:

- Crumb cup strainers and tail pieces for each bay to connect to the facility drain.
- Hot and cold water fixture that mounts to the backsplash, one for double sink, two for triple sink.
- Pure water fixture for connection of mineral free water or similar that mounts to the backsplash, one provided per Clean Up Station.
- Wall mounting brackets that attach to backsplash and wall for ease of installation.
- Leveling feet to provide adjustment up to 1½" (38mm).

#### **ARCHITECTURAL SPECIFICATION**

Material: Type 304 Stainless Steel, 16 gauge Welds: HeliArc welded and polished

Hardware: Water Fixtures by WaterSaver

Undercoating: Applied to all unexposed areas, reduces


noise and condensation build up

Finish: #4 satin finish


# 

### CLEAN UP STATION 2 SINK

S/S SINKTOP & BACKSPLASH



## **CLEAN UP STATION** 3 SINK



S/S SINKTOP & BACKSPLASH

